CALCIUM ALGINATE WOUND DRESSING

PRODUCT DESCRIPTION

Product is a disposable sterile product, which is consisted of calcium alginate fiber, it is produced by non - weaving process, and it has high absorbency and hemostatic function, which has inhibitory effect **on bacteria. The product also creates and maintains a moist wound healing environment, and reduces the secondary damage when change wound dressing.**


STRUCTURE PERFORMANCE

1) Product composition and structure

The main raw material of the products is calcium alginate fiber. Ordinary type dressing is mainly consisted of calcium alginate fiber. Island type dressing is mainly consisted of three layers of structure, the bottom film is PU/ spunlace tape, calcium alginate non-woven fabric, and PET release film:

2) Product performance description

- a、 High absorption
- b、 Promote the wound healing
- c、 Reducing scar
- d、 Hemostasis

When contact with the exudates, there are a large number of hydrophilic groups in the calcium alginate fiber, the product can absorb the water from the blood and secretion quickly and extensively. Ca^{2+} in calcium alginate and Na^+ in the secretion are exchanged, forming soft, sticky, moist alginate gel. It is not easy to lose the moisture and providing a moist healing environment for the wound surface. According to the wet environmental healing theory, maintaining a moist wound healing environment can stimulate the growth and regeneration of capillary, keeping fresh wound, accelerating the growth of granulation tissue and the formation of the epithelium can promote wound healing.

FUNDAMENTAL PRINCIPLE

- 1) Wash the wound with normal saline.
- 2) After the skin around the wound is dry, select the suitable size of the product, and the product should cover the wound completely.
- 3) If you use the island type product, tear off the release film and apply the product to the wound, then tear the back lining on the PU tape and close the PU tape to the skin. If you use ordinary product, please fix with secondary dressing.

It is suitable for all type of the wound with large exudates, keep the wound moist and absorb the exudate.

METHOD OF USE

- 1) Do not used for 3 degree burns.
- 2) Do not used for uncontrolled bleeding wounds.
- 3) Patients who are known to be allergic to calcium alginate are forbidden, cautiously use of allergic constitution.
- 4) Dry wound is not applicable.

SCOPE OF APPLICATION

- 1) This product is a disposable sterile product. It is strictly prohibited to reuse the product. Do not use if package is damaged.
 - 2) Please use it within the validity period. Read the instructions carefully before use.
- Before use of this product for refractory wound, should make a scientific analysis of wound types and causes and formulate corresponding treatment plans. Only this product can not be used as a refractory wound treatment scheme.

CONTRAIN-DICATIONS

Non-corrosive gas, cool, well ventilated, clean room, room temperature preservation.

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